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- VivaglassCEM (日本未導入)**
VivaglassCEM PL (日本未導入) is a highly translucent self-curing, radiopaque glass-ionomer cement.

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- Telio CS Inlay (日本未導入)**
Temporary light-curing filling material for deep inlay preparations with parallel walls and sealing of implant screw access holes.

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
- Cervitec Plus (日本未導入)**
The protective varnish containing chlorhexidine and thymol protects exposed root surfaces and controls bacteria.

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The screw access openings are thoroughly rinsed with water spray and dried with oil-free air. Subsequently, the screw access openings are sealed with cotton wool or foam pellets and **Telio CS Inlay** (日本未導入) . For all further treatments steps, relative isolation of the operating field, e.g. with **OptraGate**, indispensable. A retraction cord may optionally be placed.

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
5 Vivaglass CEM (日本未導入) is mixed and applied

The **Vivaglass CEM** (日本未導入) powder and liquid are mixed in a 1:1 ratio. For the cementation of a bridge, the amount of material is increased depending on the number of abutment teeth involved.

The luting material is applied to the restoration with a spatula or brush.

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Once the cement has completely set, excess is removed e.g. with an implant scaler. The setting time is 4-6 minutes.

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A thin layer of **Cervitec Plus** (日本未導入) is applied where it is needed with the help of a Vivadent applicator or a brush. The varnish sets by itself or with the application of a stream of air.